


LIFT SLINGS

SIT-TO-STAND BUTTOCK STRAP

Features

- Designed for use with the LF2020 STS lift
- Polyester with nylon stitching and webbing
- Maximum Weight Capacity
400 lb - EVENLY DISTRIBUTED
- Six (6) month limited warranty



	Size	Torso Girth		total Width		body Length		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	
DSLSA9	Buttock Strap	34.25	87	7.5	19.05	71.25	180.98	400	181.4	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF070007/RevC23