

INFECTION CONTROL

EYE GOGGLES

Features

- Splash guard applications
- Indirect ventilation
- Protects against flying fragments, objects, large chips, particles, sand or dirt
- PVC frame
- Adjustable elastic webbing headband contains latex
- Meets requirements of ANSI Z87.1

Caution: This product contains natural rubber latex, which may cause allergic reactions.



		UOM
9675	One Size Fits All	24/bx



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2013, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.