



Set n' Go Height Adjustable Rollator

The Lumex Set n' Go Height Adjustable Rollator offers the convenience and flexibility of a two-in-one product all rolled into one. This accommodates for a broad range of user heights which in turn offers more comfort for the user when seated.

- Easy height adjustment on the frame allows seat height to be positioned from 18" to 22" in 1" increments
- Comfortable, wide, seamless padded seat with zippered pouch under seat for security of personal items
- Also includes black pouch for storage of larger items
- Adjustable handle height with knobs to secure in place
- Easy to operate locking and "quick stop" loop brakes
- 6" wheel for both indoor and outdoor use
- 300 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Designed to meet the requirements of HCPCS: E0143 and E0156
- Limited lifetime warranty

Sure-Gait™ II Three Wheeled Steel Rollator Model 609101A

As a steel frame, three wheel rollator, the Sure-Gait II provides maximum maneuverability and comfort for those who do not require the extra stability of a four wheel rollator.

- Chrome-plated steel frame
- Adjustable handle height
- Easy one handed folding design
- Locking loop brakes
- Solid, puncture resistant 8" wheels for indoor and outdoor use
- Carry pouch included
- 300 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Designed to meet the requirements of HCPCS code E0143



GF 3-Wheel Cruiser Model 609201

The GF 3-Wheel Cruiser is designed to provide convenient mobility for those individuals who do not require a seat or the additional stability of a four-wheel rollator. This new light weight aluminum walking aid allows the user to have more freedom in movement to perform daily activities.

- Lightweight aluminum frame - weighs only 10.5 lbs
- Ergonomic hand grips with locking loop brakes
- Adjustable handle height for broad range of users
- 8" black wheels for indoor and outdoor use
- Comes complete with zippered carry pouch
- 250 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Available in 3 new attractive colors: Pewter Metallic, Metallic Burgundy and Bondi Blue
- Designed to meet the requirements of HCPCS code E0143

Deluxe Four-Wheeled Rollator Model RJ4200A

A versatile walker with four oversized 8" wheels that provide a soft ride over most types of terrain.

- Heavy-duty powder coated steel frame
- Ergonomic hand grips with locking loop brakes
- Standard with shopping basket and removable tray
- Foldable for easy storage and transport
- 500 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Designed to meet the requirements of HCPCS codes E0149 and E0156
- Includes laminated color operating instructions



GRAHAM-FIELD

For more information on the Lumex Walkabout Rollators, visit www.grahamfield.com or call 770.368.4700 © 2012, GF Health Products, Inc.

GF1000097 RevC12

Get What
You Need... & A Whole Lot More From Lumex®,
 The Most Trusted Name In Walking Aids.

LUMEX® For The Quality Of Life



GRAHAM-FIELD



The Lumex Walkabout™ Four-Wheel Rollators provide the ideal combination of comfort, appeal, and affordability. The Walkabout series offers both high quality and a variety of options to meet the needs of a wide spectrum of patients. No two patients are alike, so let Lumex, with our complete line of rollators, help put their lives back in motion.



As part of our Quality initiative, our walkabout series meets ISO standards.



Model RJ4400

Model RJ4405

Model RJ4300

Model RJ4805

Model RJ4302

Model RJ4301

Lumex® Rollator Features

Designed with an aluminum frame built for durability, these lightweight rollators allow for easy maneuverability and fold quickly and effortlessly for transport. Each rollator features a wide padded seat and tool-free removable padded backrest for enhanced user comfort. All models have ergonomic hand grips that provide a better grip as well as locking loop brakes that are easy to operate and locking rear wheels for stable positioning when sitting.

The adjustable handle height accommodates a broad range of user heights with hemi models available for shorter individuals. Each comes with a wire basket that can be placed underneath the seat or, on some models, in front for carrying personal items. The Walkabout series is designed to meet the requirements of HCPCS codes E0143 and E0156 (E0149 for Bariatric models) and offers a Limited Lifetime Warranty as stated on the laminated color operating instructions included.

Walkabout™ Imperial Bariatric Models RJ4400, RJ4405

The Walkabout Imperial rollator is designed for the bariatric user, featuring a 20"-wide padded seat and offers a 500 lb maximum weight capacity EVENLY DISTRIBUTED. It also features heavy-duty 8" wheels that are ideal for both indoor and outdoor rugged terrain. Offered in both standard and ConTour options.

Walkabout™ Deluxe ConTour Model RJ4805

The Walkabout Deluxe ConTour rollators offer the same benefits as the Walkabout Lite rollator, while featuring a contoured padded backrest and heavy-duty 8" wheels that are ideal for both indoor and outdoor rugged terrain. Also available is an add on adjustable backrest (RJ4805-ABAR) that allows the user to adjust the backrest for ideal positioning.

Walkabout™ Lite Model RJ4300

Our most popular model weighs only 14.5 lbs and is offered in various attractive colors: teal green, burgundy, pink, royal blue, aqua, lavender, champagne, and black. Each has 6" wheels for indoor or outdoor use and has a 300 lb. maximum weight capacity EVENLY DISTRIBUTED. The walkabout lite folds effortlessly with folding strap for convenient transport and storage.

Walkabout Wide Model RJ4318

Designed for the user who is "in between sizes" of a standard or bariatric rollator, this Walkabout Wide offers more room and added comfort. The wider frame features an 18.5" wide opening in between the handgrips, a 15.5" wide padded seat and a 350 lb maximum weight capacity EVENLY DISTRIBUTED.



RJ4318 offers 18.5" between handgrips



Easy grip knob makes for simple backbar removal for storage and transport.

Available colors
*Color offerings vary for different models.



New Rear Reflectors increase visibility at night!

Walkabout™ Lite Hemi Models RJ4302, RJ44002

This Hemi rollator offers the same benefits as the Walkabout Lite rollator yet is designed with a lower seat to floor height to accommodate shorter individuals. 6" wheels are available on this model for both indoor and outdoor use. An Imperial Bariatric Hemi is also available (RJ4402R) to support a 500 lb. maximum weight capacity EVENLY DISTRIBUTED.

Walkabout™ Junior Model RJ4301

This lightweight Walkabout Junior rollator is ideal for the petite or junior user. It is designed with a smaller frame, weighs only 12.5 lbs and has 6" wheels for indoor and outdoor use.